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GENERAL CONSIDERATIONS:

- The D29 implant system is a mixed system using not only the properties of a screw implant, but also the possibility of stabilisation of the implant with a plate screwed above the bone when the implantation is made on a poor quality or too-thin bone.
- To ensure it is tolerated in the bone and in the tissue, the D29 SIT system is manufactured in titanium according to the ISO standards for implantable-grade titanium.
- The D29 SIT shall always be implanted by an experienced professional trained to SIT implant techniques.

WARNING:

- The D29 implant system is composed of specific additional parts that cannot be replaced by any other brand than SIT. If you use any prosthetic material and/or fastenings not validated by SIT, the manufacturer warranty will *de facto* expire.
- In no case should the implantable parts of D29 systems be reused (single use).
- The stabilisation plate fixation in a cortical bone, determined at the practitioner's discretion, is not intended in areas where the absence of cortical bone or its poor quality contraindicates it. This fixation will be exclusively performed with SIT screws. Its integration depends on the quality of the plate positioning and the fixation screwed in the cortical bone.

PRECAUTIONS:

- To install implants from the D29 line, it is essential that the practitioner took a specific training exclusively proposed by SIT in order to master the various items to use taking into account the patient's anatomical constraints.
- The 2-mm-diameter screws provided with the D29 plate kit are "self-drilling" screws and are installed without any prior drilling. Nevertheless, in presence of exceptionally dense cortical bone, it could be necessary to drill beforehand the bone with the adapted drill.

INDICATIONS:

- The D29 SIT system is a mixed system composed of a screw implant and a stabilisation plate that can be used if necessary if the bone characteristics do not allow any primary stabilisation appropriate for the immediate loading of a dental prosthesis or, at the practitioner's discretion, the increased stability of the screw implant.
- As the D29 SIT stabilisation plate is removable, the practitioner, at his/her discretion, can remove it, at a distance from the implantation.

CONTRAINDICATIONS:

- Product: In no case should the D29 SIT system be installed in areas where its screw implant part cannot be anchored, either because of insufficient height and/or bone density.
- Patient: Local or systemic chronic or acute infection / Poor oral hygiene / Significant bruxism.
- Neuromuscular and/or mental disorders that may lead to the patient's non-respect of the precautions recommended by the practitioner.
- A very special attention should be turned to patients that may have been irradiated.

SIDE EFFECTS:

- All the side effects that can occur after surgical breaking and the installation of one or several implants (oedema – haemorrhage – inflammation – paresthesia...)

PATIENT INFORMATION:

- The practitioner installing the D29 SIT should inform the patient that the implantation safety and/or durability depends on his/her activity and on the respect of the recommendations he/she will have received.
- After the implantation, the patient should indicate to the surgeon any sensation or modification appearing at the level of the D29 SIT implant.
- In case of treatment or investigation (MRI, for example), the patient shall inform the practitioner about the presence of the D29 SIT.

ANCILLARY MATERIAL:

- Only use the instruments and various ancillary materials provided in the SIT D29 kit to install the implants of the D29 system.

PACKAGING, STORAGE:

- The implantable D29 SIT system is provided sterile in a unit packaging.
- Every packaging has, in addition to the product references, self-adhesive labels stating all the manufacturing specifications (batch number – sterilisation – expiration date, etc.), which will be noted in the patient file.
- In no case should any implant taken out from its packaging, thus not sterile anymore, be reused.

DECONTAMINATION, CLEANING:

- When the implants are not sterile anymore, once taken out from the packaging, they should be decontaminated, and then cleaned.
- Particularly before sterilising the ancillary materials, a decontamination stage shall be performed according to the current standards (soaking into an adapted liquid, followed by rinsing into plenty of water colder than 30° C.
- Cleaning is performed with an appropriate liquid making sure the products are not being damaged (ancillary or implants). The cleaning stage can be performed manually, in a traditional washing machine or with ultrasounds.
- In any case, the screw implants or plates shall never be sterilised again when they are taken out from their initial packaging.

PACKAGING BEFORE STERILISATION:

- Immediately performed after unpacking, the cleaning and drying stages will be performed according to the current standards.

STERILISATION:

- Performed after the previous stages, it will strictly respect the current directives (for example: water steam sterilisation between 134°C and 138°C for at least 18 minutes).

USE ENVIRONMENT OF THE IMPLANTS DURING THE PROCEDURE:

- Aseptic environment according to the current standards.

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